



TSPAN6 mouse mAb

Catalog No	YP-mAb-18132
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TSPAN6 TM4SF6 UNQ767/PRO1560
Immunogen	Synthesized peptide derived from human TSPAN6
Specificity	This antibody detects endogenous levels of TSPAN6 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	A15 homolog;T245;T245 protein;tetraspanin 6;Tetraspanin TM4 D;TM4SF6;TSPAN 6;TSPAN6;TSPAN6 TM4SF6 UNQ767/PRO1560;Tetraspanin-6 (Tspan-6) (A15 homolog) (Putative NF-kappa-B-activating protein 321) (T245 protein) (Tetraspanin TM4-D) (Transmembrane 4 superfamily member 6)
Observed Band	25-30 kDa
Calculated Molecular Weight	245 aa, 28 kDa
Cell Pathway	Membrane; Multi-pass membrane protein.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images