



FBXO28 mouse mAb

Catalog No	YP-mAb-18138
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FBXO28 KIAA0483
Immunogen	Synthesized peptide derived from human FBXO28
Specificity	This antibody detects endogenous levels of FBXO28 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F box only protein 28;F box protein 28;Fbx28;FBXO28;KIAA0483;FBXO28 KIAA0483;F-box only protein 28
Observed Band	42 kDa
Calculated Molecular Weight	368 aa, 41 kDa
Cell Pathway	Chromosome, centromere, kinetochore .
Tissue Specificity	
Function	Probably recognizes and binds to some phosphorylated proteins and promotes their ubiquitination and degradation.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images