



LRRC1 mouse mAb

Catalog No	YP-mAb-18139
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LRRC1 LANO
Immunogen	Synthesized peptide derived from human LRRC1
Specificity	This antibody detects endogenous levels of LRRC1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Leucine-rich repeat-containing protein 1;LAP and no PDZ protein;LANO adapter protein;LANO;LRRC1;LRRC1 LANO;Leucine-rich repeat-containing protein 1 (LANO adapter protein) (LAP and no PDZ protein)
Observed Band	59 kDa
Calculated Molecular Weight	59 kDa
Cell Pathway	Cytoplasm . Membrane ; Peripheral membrane protein . Localized at the basolateral side of epithelial cells.
Tissue Specificity	Expressed strongly in testis and placenta, followed by heart, lung, kidney, thyroid, trachea, colon, prostate and pancreas.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images