



ATGL mouse mAb

Catalog No	YP-mAb-18140
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PNPLA2 ATGL FP17548
Immunogen	Synthesized peptide derived from human ATGL
Specificity	This antibody detects endogenous levels of ATGL at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Adipose triglyceride lipase;PNPLA2;Calcium-independent phospholipase A2-zeta;Desnutrin;EC:3.1.1.3;PNPLA2 ATGL FP17548;Patatin-like phospholipase domain-containing protein 2 (EC 3.1.1.3) (Adipose triglyceride lipase) (Calcium-independent phospholipase A2) (Desnutrin) (IPLA2-zeta) (Pigment epithelium-derived factor)
Observed Band	45-55 kDa
Calculated Molecular Weight	55 kDa
Cell Pathway	Lipid droplet . Cell membrane ; Multi-pass membrane protein . Cytoplasm .
Tissue Specificity	Highest expression in adipose tissue. Also detected in heart, skeletal muscle, and portions of the gastrointestinal tract. Detected in normal retina and retinoblastoma cells. Detected in retinal pigment epithelium and, at lower intensity,
Function	Catalyzes the initial step in triglyceride hydrolysis in adipocyte and non-adipocyte lipid droplets . Exhibits a strong preference for the hydrolysis of long-chain fatty acid esters at the sn-2 position of the glycerol backbone and acts coordinately with LIPE/HLS and DGAT2 within the lipolytic cascade (By similarity). Also possesses acylglycerol transacylase and phospholipase A2 activities . Transfers fatty acid from triglyceride to retinol, hydrolyzes



retinylesters, and generates 1,3-diacylglycerol from triglycerides . Regulates adiposome size and may be involved in the degradation of adiposomes . May play an important role in energy homeostasis (By similarity). May play a role in the response of the organism to starvation, enhancing hydrolysis of triglycerides and providing free fatty acids to other tissues to be oxidized in situations of energy depletion (By similarity).

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images