



PRPS2 mouse mAb

Catalog No	YP-mAb-18148
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PRPS2
Immunogen	Synthesized peptide derived from human PRPS2
Specificity	This antibody detects endogenous levels of PRPS2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:2.7.6.1;PPRibP;PRS II;PRSI;PRS-II;PRPS2;Ribose-phosphate pyrophosphokinase 2 (EC 2.7.6.1) (PPRibP) (Phosphoribosyl pyrophosphate synthase II) (PRS-II)
Observed Band	3-34 kDa
Calculated Molecular Weight	35kD
Cell Pathway	
Tissue Specificity	
Function	Catalyzes the synthesis of phosphoribosylpyrophosphate (PRPP) that is essential for nucleotide synthesis.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images