



PLAC8 mouse mAb

Catalog No	YP-mAb-18150
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PLAC8 BM-004
Immunogen	Synthesized peptide derived from human PLAC8
Specificity	This antibody detects endogenous levels of PLAC8 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BM-004;C15;onzin;placenta 8;Placenta-specific gene 8 protein;PLAC8;PLAC8 BM-004;Placenta-specific gene 8 protein (Protein C15)
Observed Band	12 kDa
Calculated Molecular Weight	12 kDa
Cell Pathway	
Tissue Specificity	Expressed at high levels in plasmacytoid dendritic cells. High expression in spleen, lymph nodes, peripheral blood leukocytes, and bone marrow, with lower expression in thymus, appendix, and fetal liver.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images