



CHCHD5 mouse mAb

Catalog No	YP-mAb-18158
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CHCHD5 C2orf9
Immunogen	Synthesized peptide derived from human CHCHD5
Specificity	This antibody detects endogenous levels of CHCHD5 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C2orf9;Coiled-coil-helix-coiled-coil-helix domain-containing protein 5;CHCHD5;CHCHD5 C2orf9
Observed Band	12 kDa
Calculated Molecular Weight	110 aa, 12 kDa
Cell Pathway	Mitochondrion intermembrane space .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images