



ARL13B mouse mAb

Catalog No	YP-mAb-18160
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ARL13B ARL2L1
Immunogen	Synthesized peptide derived from human ARL13B
Specificity	This antibody detects endogenous levels of ARL13B at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ADP-ribosylation factor-like protein 13B;ADP-ribosylation factor-like protein 2-like 1;ARL2 like protein 1;ARL2L1;ARL2-like protein 1;ARL13B;ARL13B ARL2L1;ADP-ribosylation factor-like protein 13B (ADP-ribosylation factor-like protein 2-like 1) (ARL2-like protein 1)
Observed Band	40-48 kDa, 60 kDa
Calculated Molecular Weight	48 kDa
Cell Pathway	Cell projection, cilium membrane ; Lipid-anchor . Cell projection, cilium . Associates to the cilium membrane via palmitoylation. Localizes to proximal ciliary membranes, to an inversin-like subciliary membrane compartment, excluding the transition zone.
Tissue Specificity	Expressed in the developing brain.
Function	Cilium-specific protein required to control the microtubule-based, ciliary axoneme structure. May act by maintaining the association between IFT subcomplexes A and B. Binds GTP but is not able to hydrolyze it; the GTPase activity remains unclear. Required to pattern the neural tube. Involved in cerebral cortex development: required for the initial formation of a polarized radial glial scaffold, the first step in the construction of the cerebral cortex, by regulating ciliary signaling. Regulates the migration and placement of postmitotic



interneurons in the developing cerebral cortex. May regulate endocytic recycling traffic; however, additional evidence is required to confirm these data.

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images