



COMMD8 mouse mAb

Catalog No	YP-mAb-18171
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	COMMD8 MDS022
Immunogen	Synthesized peptide derived from human COMMD8
Specificity	This antibody detects endogenous levels of COMMD8 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	COMM domain-containing protein 8;COMMD 8;MDS022;COMMD8;COMMD8 MDS022
Observed Band	21 kDa
Calculated Molecular Weight	183 aa, 21 kDa
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	Widely expressed with highest expression in thyroid.
Function	May modulate activity of cullin-RING E3 ubiquitin ligase (CRL) complexes . May down-regulate activation of NF-kappa-B .
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images