



ANKRD34B mouse mAb

Catalog No	YP-mAb-18180
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ANKRD34B
Immunogen	Synthesized peptide derived from human ANKRD34B
Specificity	This antibody detects endogenous levels of ANKRD34B at Human, Mous
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ANKRD34B;ankyrin repeat domain 34B;DP58;Ankyrin repeat domain-containing protein 34B
Observed Band	55-66 kDa
Calculated Molecular Weight	514 aa, 56 kDa
Cell Pathway	Cytoplasm. Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images