



TWIST2 mouse mAb

Catalog No	YP-mAb-18184
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TWIST2 BHLHA39 DERMO1
Immunogen	Synthesized peptide derived from human TWIST2
Specificity	This antibody detects endogenous levels of TWIST2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	bHLHa39;Dermis expressed protein 1;Dermo 1;DERMO1;twist homolog 2 (Drosophila);Twist related protein 2;TWIST2;TWIST2 BHLHA39 DERMO1;Twist-related protein 2 (Class A basic helix-loop-helix protein 39) (bHLHa39) (Dermis-expressed protein 1) (Dermo-1)
Observed Band	21 kDa, 60 kDa
Calculated Molecular Weight	160 aa, 18 kDa
Cell Pathway	Nucleus . Cytoplasm . Mainly nuclear during embryonic development. Cytoplasmic in adult tissues.
Tissue Specificity	In the embryo, highly expressed in chondrogenic cells. In embryonic skin, expressed in the undifferentiated mesenchymal layer beneath the epidermis which later develops into the dermis. Expressed in early myeloid cells but not in lymphoid cells in the liver. Expression also detected in the secretory ependymal epithelium of the choroid plexus primordium. In the adult, expressed in secreting glandular tissues and
Function	Binds to the E-box consensus sequence 5'-CANNTG-3' as a heterodimer and inhibits transcriptional activation by MYOD1, MYOG, MEF2A and MEF2C. Also represses expression of pro-inflammatory cytokines such as TNFA and IL1B. Involved in postnatal glycogen storage and energy metabolism (By similarity).



Inhibits the premature or ectopic differentiation of preosteoblast cells during osteogenesis, possibly by changing the internal signal transduction response of osteoblasts to external growth factors.

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images