



SNX25 mouse mAb

Catalog No	YP-mAb-18189
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SNX25 MSTP043
Immunogen	Synthesized peptide derived from human SNX25
Specificity	This antibody detects endogenous levels of SNX25 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Sorting nexin-25;sorting nexin 25;SBBI31;MSTP043;SNX25;SNX25 MSTP043
Observed Band	98-105 kDa, 65 kDa
Calculated Molecular Weight	98 kDa
Cell Pathway	May be involved in several stages of intracellular trafficking.
Tissue Specificity	
Function	May be involved in several stages of intracellular trafficking.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images