



CTHRC1 mouse mAb

Catalog No	YP-mAb-18199
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CTHRC1 UNQ762/PRO1550
Immunogen	Synthesized peptide derived from human CTHRC1
Specificity	This antibody detects endogenous levels of CTHRC1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	collagen triple helix repeat containing 1;Collagen triple helix repeat-containing protein 1;Protein NMTC1;UNQ762/PRO1550;CTHRC1;CTHRC1 UNQ762/PRO1550;Collagen triple helix repeat-containing protein 1 (Protein NMTC1)
Observed Band	25-30 kDa
Calculated Molecular Weight	243 aa, 26 kDa
Cell Pathway	Secreted, extracellular space, extracellular matrix .
Tissue Specificity	Isoform 1 is expressed in calcified atherosclerotic plaque and chondrocyte-like cells.
Function	May act as a negative regulator of collagen matrix deposition.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images