



C1orf93 mouse mAb

Catalog No	YP-mAb-18207
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FAM213B C1orf93
Immunogen	Synthesized peptide derived from human C1orf93
Specificity	This antibody detects endogenous levels of C1orf93 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C1orf93;FAM213B;RP3 395M2.1;PRXL2B;FAM213B C1orf93;Prostamide/prostaglandin F synthase (Prostamide/PG F synthase) (Prostamide/PGF synthase) (EC 1.11.1.20) (Protein FAM213B)
Observed Band	21 kDa
Calculated Molecular Weight	22kD
Cell Pathway	Cytoplasm, cytosol .
Tissue Specificity	
Function	Catalyzes the reduction of prostaglandin-ethanolamide H(2) (prostamide H(2)) to prostamide F(2alpha) with NADPH as proton donor. Also able to reduce prostaglandin H(2) to prostaglandin F(2alpha) (By similarity).
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images