



SNRPB2 mouse mAb

Catalog No	YP-mAb-18215
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SNRPB2
Immunogen	Synthesized peptide derived from human SNRPB2
Specificity	This antibody detects endogenous levels of SNRPB2 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	U2 snRNP B";U2 snRNP B;U2 small nuclear ribonucleoprotein B";SNRPB2;U2 small nuclear ribonucleoprotein B" (U2 snRNP B")
Observed Band	25-31 kDa
Calculated Molecular Weight	225 aa, 26, 5 kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	Involved in pre-mRNA splicing as component of the spliceosome . Associated with sn-RNP U2, where it contributes to the binding of stem loop IV of U2 snRNA .
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images