



VAMP3/Cellubrevin mouse mAb

Catalog No	YP-mAb-18220
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	VAMP3 SYB3
Immunogen	Synthesized peptide derived from human VAMP3/Cellubrevin
Specificity	This antibody detects endogenous levels of VAMP3/Cellubrevin at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	VAMP3;CEB;Cellubrevin;SYB3;Synaptobrevin-3;VAMP3 SYB3;Vesicle-associated membrane protein 3 (VAMP-3) (Cellubrevin) (CEB) (Synaptobrevin-3)
Observed Band	11-17 kDa
Calculated Molecular Weight	11 kDa
Cell Pathway	Membrane ; Single-pass type IV membrane protein . Cell junction, synapse, synaptosome .
Tissue Specificity	
Function	SNARE involved in vesicular transport from the late endosomes to the trans-Golgi network.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images