



MIPEP mouse mAb

Catalog No	YP-mAb-18223
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MIPEP MIP
Immunogen	Synthesized peptide derived from human MIPEP
Specificity	This antibody detects endogenous levels of MIPEP at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HMIP;MIP;MIPEP;MIPEP MIP;Mitochondrial intermediate peptidase (MIP) (EC 3.4.24.59)
Observed Band	70-80 kDa
Calculated Molecular Weight	81 kDa
Cell Pathway	Mitochondrion matrix.
Tissue Specificity	
Function	Cleaves proteins, imported into the mitochondrion, to their mature size.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images