



SDHA mouse mAb

Catalog No	YP-mAb-18224
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SDHA
Immunogen	A synthesized peptide derived from human SDHA
Specificity	Affinity Chromatography
Formulation	Unconjugated
Source	Unmodified
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:1.3.5.1;FP;MC2DN1;PGL5;PPGL5;SDHA;SDH2;SDHF;Succinate dehydrogenase [ubiquinone] flavoprotein subunit;mitochondrial;Flavoprotein subunit of complex II
Observed Band	Calculated MW: 73 kDa; Observed MW: 73 kDa
Calculated Molecular Weight	70kD
Cell Pathway	
Tissue Specificity	
Function	
Background	Flavoprotein (FP) subunit of succinate dehydrogenase (SDH) that is involved in complex II of the mitochondrial electron transport chain and is responsible for transferring electrons from succinate to ubiquinone.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images