



LC3B mouse mAb

Catalog No	YP-mAb-18227
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MAP1LC3B
Immunogen	Recombinant Protein of MAP LC3 β
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	$\geq 90\%$
Storage Stability	-20°C/1 year
Synonyms	MAP1LC3B;ATG8F;Autophagy-related ubiquitin-like modifier LC3 B;LC3B;MAP1 light chain 3-like protein 2
Observed Band	Calculated MW: 15 kDa; Observed MW: 14,16 kDa
Calculated Molecular Weight	Calculated MW: 15 kDa; Observed MW: 14,16 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

