



APOBEC3B mouse mAb

Catalog No	YP-mAb-18233
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	APOBEC3B
Immunogen	Synthesized peptide derived from human APOBEC3B
Specificity	This antibody detects endogenous levels of APOBEC3B at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	A3B;DNA dC->dU-editing enzyme APOBEC-3B;EC:3.5.4.38;Phorbolin 1 related protein;Phorbolin 2/3;APOBEC3B;Probable DNA dC->dU-editing enzyme APOBEC-3B (EC 3.5.4.-) (Phorbolin-1-related protein) (Phorbolin-2/3)
Observed Band	46 kDa
Calculated Molecular Weight	46 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Expressed at high and moderate levels in peripheral blood leukocytes, spleen, testes, heart, thymus, prostate and ovary. Also expressed at low levels in several other tissues.
Function	DNA deaminase (cytidine deaminase) which acts as an inhibitor of retrovirus replication and retrotransposon mobility via deaminase-dependent and - independent mechanisms. After the penetration of retroviral nucleocapsids into target cells of infection and the initiation of reverse transcription, it can induce the conversion of cytosine to uracil in the minus-sense single-strand viral DNA, leading to G-to-A hypermutations in the subsequent plus-strand viral DNA. The resultant detrimental levels of mutations in the proviral genome, along with a deamination-independent mechanism that works prior to the proviral integration, together exert efficient antiretroviral effects in infected target cells. Selectively



targets single-stranded DNA and does not deaminate double-stranded DNA or single- or double-stranded RNA. Exhibits antiviral activity against simian immunodeficiency virus (SIV), hepatic

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images