



SAMD3 mouse mAb

Catalog No	YP-mAb-18239
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SAMD3
Immunogen	Synthesized peptide derived from human SAMD3
Specificity	This antibody detects endogenous levels of SAMD3 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SAMD3;Sterile alpha motif domain-containing protein 3 (SAM domain-containing protein 3)
Observed Band	52 kDa
Calculated Molecular Weight	221aa,25 kDa; 520aa,61 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images