



GOLPH3 mouse mAb

Catalog No	YP-mAb-18240
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GOLPH3 GPP34
Immunogen	Synthesized peptide derived from human GOLPH3
Specificity	This antibody detects endogenous levels of GOLPH3 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MIDAS;GPP34;GOPP1;GOLPH 3;Golgi phosphoprotein 3;GOLPH3;GOLPH3 GPP34;Golgi phosphoprotein 3 (Coat protein GPP34) (Mitochondrial DNA absence factor) (MIDAS)
Observed Band	34 kDa
Calculated Molecular Weight	298 aa, 34 kDa
Cell Pathway	Golgi apparatus, Golgi stack membrane; Peripheral membrane protein; Cytoplasmic side. Golgi apparatus, trans-Golgi network membrane; Peripheral membrane protein; Cytoplasmic side. Mitochondrion intermembrane space. Cell membrane . Endosome . Phosphatidylinositol 4-phosphate-binding and oligomerization participate in the recruitment onto Golgi membranes. .
Tissue Specificity	Detected in muscle fibers of patients with mitochondrial diseases; not detected in normal muscle fibers.
Function	Phosphatidylinositol-4-phosphate-binding protein that links Golgi membranes to the cytoskeleton and may participate in the tensile force required for vesicle budding from the Golgi. Thereby, may play a role in Golgi membrane trafficking and could indirectly give its flattened shape to the Golgi apparatus. May also bind to the coatomer to regulate Golgi membrane trafficking. May play a role in anterograde transport from the Golgi to the plasma membrane and regulate



secretion. Has also been involved in the control of the localization of Golgi enzymes through interaction with their cytoplasmic part. May play an indirect role in cell migration. Has also been involved in the modulation of mTOR signaling. May also be involved in the regulation of mitochondrial lipids biosynthesis.

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images