



GOLGA3 mouse mAb

Catalog No	YP-mAb-18255
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GOLGA3
Immunogen	Synthesized peptide derived from human GOLGA3
Specificity	This antibody detects endogenous levels of GOLGA3 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GCP170;Golgi complex-associated protein of 170 kDa;Golgin 160;Golgin-160; MEA 2;GOLGA3;Golgin subfamily A member 3 (Golgi complex-associated protein of 170 kDa) (GCP170) (Golgin-160)
Observed Band	167 kDa
Calculated Molecular Weight	1498 aa, 167 kDa
Cell Pathway	Cytoplasm. Golgi apparatus, Golgi stack membrane; Peripheral membrane protein.
Tissue Specificity	Expressed in all tissues tested. Expressed in liver, testis, lung, heart, salivary gland and kidney.
Function	Golgi auto-antigen; probably involved in maintaining Golgi structure.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images