



SNRNP200 mouse mAb

Catalog No	YP-mAb-18257
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SNRNP200 ASCC3L1 HELIC2 KIAA0788
Immunogen	Synthesized peptide derived from human SNRNP200
Specificity	This antibody detects endogenous levels of SNRNP200 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Activating signal cointegrator 1 complex subunit 3-like 1;ASCC3L1;BRR2; EC:3.6.4.13;HELIC2;SNRNP200;SNRNP200 ASCC3L1 HELIC2 KIAA0788;U5 small nuclear ribonucleoprotein 200 kDa helicase (EC 3.6.4.13) (Activating signal cointegrator 1 complex subunit 3-like 1) (BRR2 homolog) (U5 snRNP-specific 200 kDa protein) (U5-200KD)
Observed Band	200-245 kDa
Calculated Molecular Weight	494 aa, 57 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Widely expressed.
Function	Plays role in pre-mRNA splicing as core component of precatalytic, catalytic and postcatalytic spliceosomal complexes . Involved in spliceosome assembly, activation and disassembly. Mediates changes in the dynamic network of RNA-RNA interactions in the spliceosome. Catalyzes the ATP-dependent unwinding of U4/U6 RNA duplexes, an essential step in the assembly of a catalytically active spliceosome.
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images