



NR2F1 mouse mAb

Catalog No	YP-mAb-18260
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NR2F1 EAR3 ERBAL3 TFCOUP1
Immunogen	Synthesized peptide derived from human NR2F1
Specificity	This antibody detects endogenous levels of NR2F1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	COUP TF I;COUP TF1;COUP TFI;COUP transcription factor 1;COUP transcription factor I;NR2F1;NR2F1 EAR3 ERBAL3 TFCOUP1;COUP transcription factor 1 (COUP-TF1) (COUP transcription factor I) (COUP-TF I) (Nuclear receptor subfamily 2 group F member 1) (V-erbA-related protein 3) (EAR-3)
Observed Band	46 kDa
Calculated Molecular Weight	50 kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	Coup (chicken ovalbumin upstream promoter) transcription factor binds to the ovalbumin promoter and, in conjunction with another protein (S300-II) stimulates initiation of transcription. Binds to both direct repeats and palindromes of the 5'-AGGTCA-3' motif. Represses transcriptional activity of LHCG.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images