



DIMT1L mouse mAb

Catalog No	YP-mAb-18263
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DIMT1 DIMT1L HUSSY-05
Immunogen	Synthesized peptide derived from human DIMT1L
Specificity	This antibody detects endogenous levels of DIMT1L at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	18S rRNA dimethylase;DIM1;DIMT1;DIMT1L;HSA9761;HUSSY5;DIMT1 DIMT1L HUSSY-05;Probable dimethyladenosine transferase (EC 2.1.1.183) (DIM1 dimethyladenosine transferase 1 homolog) (DIM1 dimethyladenosine transferase 1-like) (Probable 18S rRNA (adenine(1779)-N(6)/adenine(1780)-N(
Observed Band	35 kDa
Calculated Molecular Weight	35 kDa
Cell Pathway	Nucleus, nucleoplasm . Nucleus, nucleolus .
Tissue Specificity	
Function	Specifically dimethylates two adjacent adenosines in the loop of a conserved hairpin near the 3'-end of 18S rRNA in the 40S particle . Involved in the pre-rRNA processing steps leading to small-subunit rRNA production independently of its RNA-modifying catalytic activity .
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images