



RPS14 mouse mAb

Catalog No	YP-mAb-18272
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RPS14 PRO2640
Immunogen	Synthesized peptide derived from human RPS14
Specificity	This antibody detects endogenous levels of RPS14 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Small ribosomal subunit protein uS11;RPS 14;ribosomal protein S14;PRO2640; EMTB;RPS14;RPS14 PRO2640;40S ribosomal protein S14
Observed Band	16 kDa
Calculated Molecular Weight	16 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images