



DDX41 mouse mAb

Catalog No	YP-mAb-18291
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DDX41 ABS
Immunogen	Synthesized peptide derived from human DDX41
Specificity	This antibody detects endogenous levels of DDX41 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ABS;DEAD box protein 41;EC:3.6.4.13;Probable ATP-dependent RNA helicase DDX41;DDX41;DDX41 ABS;Probable ATP-dependent RNA helicase DDX41 (EC 3.6.4.13) (DEAD box protein 41) (DEAD box protein abstract homolog)
Observed Band	70 kDa
Calculated Molecular Weight	70 kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	Probable ATP-dependent RNA helicase. Is required during post-transcriptional gene expression. May be involved in pre-mRNA splicing.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images