



L3MBTL4 Mouse mAb

Catalog No	YP-mAb-18452
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	L3MBTL4
Immunogen	Synthesized peptide derived from human L3MBTL4
Specificity	This antibody detects endogenous levels of L3MBTL4 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	H I(3)mbt like protein 4;HsT131;l(3)mbt like 4 (Drosophila);L(3)mbt like protein 4; L3MBTL4;Lethal(3)malignant brain tumor-like protein 4 (H-I(3)mbt-like protein 4) (L(3)mbt-like protein 4)
Observed Band	75-8 kDa
Calculated Molecular Weight	69kD
Cell Pathway	Nucleus .
Tissue Specificity	
Function	Putative Polycomb group (PcG) protein. PcG proteins maintain the transcriptionally repressive state of genes, probably via a modification of chromatin, rendering it heritably changed in its expressibility (By similarity).
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
