



DDX59 Mouse mAb

Catalog No	YP-mAb-18488
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human DDX59 (NP_001026895.2).
Specificity	
Formulation	
Source	
Purification	Affinity purification
Dilution	WB 1:500 - 1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DEAD box protein 59;EC:3.6.4.13;Probable ATP-dependent RNA helicase DDX59;Zinc finger HIT domain-containing protein 5;ZNHIT5;OFD5;DDX59
Observed Band	75 kDa
Calculated Molecular Weight	619 aa, 69 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Predicted to enable RNA binding activity and RNA helicase activity. Predicted to be located in cytoplasm and nucleus. Predicted to be integral component of membrane. Implicated in orofacioidigital syndrome V.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
