



UBE2J2 Monoclonal Antibody

Catalog No	YP-mAb-18808
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	UBE2J2 NCUBE2
Immunogen	Synthesized peptide derived from human UBE2J2
Specificity	This antibody detects endogenous levels of UBE2J2 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UBC6;NCUBE2;NCUBE-2;PRO2121;UBE2J2;UBE2J2 NCUBE2;Ubiquitin-conjugating enzyme E2 J2 (Non-canonical ubiquitin-conjugating enzyme 2) (NCUBE-2)
Observed Band	29kDa
Calculated Molecular Weight	29kDa
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass type IV membrane protein .
Tissue Specificity	
Function	Catalyzes the covalent attachment of ubiquitin to other proteins. Seems to function in the selective degradation of misfolded membrane proteins from the endoplasmic reticulum (ERAD).
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

