



LRRC71 C1orf92 Monoclonal Antibody

Catalog No	YP-mAb-18809
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LRRC71 C1orf92
Immunogen	Synthesized peptide derived from human C1orf92
Specificity	This antibody detects endogenous levels of C1orf92 at Human, Mouse
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NET11;CGI-49;SCCPDH;LRRC71 C1orf92;Leucine-rich repeat-containing protein 71
Observed Band	45kDa
Calculated Molecular Weight	47kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

