



RNF128 Mouse mAb

Catalog No	YP-mAb-18811
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications WB	
Gene Name	RNF128 GRAIL
Immunogen	Synthesized peptide derived from human RNF128
Specificity	This antibody detects endogenous levels of RNF128 at Human, Mouse
Formulation	
Source Monoclonal, Mouse,IgG	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	E3 ubiquitin protein ligase RNF128;E3 ubiquitin-protein ligase RNF128; FLJ23516;Gene related to anergy in lymphocytes protein;GRAIL;RING finger protein 128;RN128_HUMAN;RNF128;RNF128 GRAIL;E3 ubiquitin-protein ligase RNF128 (Gene related to anergy in lymphocytes protein) (RING finger protein 128)
Observed Band	47kD
Calculated Molecular Weight	47kD(Calculated).
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
