



ACYP1 Mouse mAb

Catalog No	YP-mAb-18922
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ACYP1 ACYPE
Immunogen	Synthesized peptide derived from human ACYP1
Specificity	This antibody detects endogenous levels of ACYP1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	erythrocyte isozyme;organ-common type isozyme;Acylphosphatase 1; Acylphosphatase 1 erythrocyte (common) type;Acylphosphatase 1 erythrocyte type;Acylphosphatase;Acylphosphatase erythrocyte isozyme;Acylphosphatase organ common type isozyme;Acylphosphatase-1;Acylphosphatase1; Acylphosphate phosphohydrolase 1;ACYP 1;ACYP1;ACYP1_HUMAN;ACYPE; Erythrocyte acylphosphatase 1;ACYP1 ACYPE;Acylphosphatase-1 (Acylphosphatase;erythrocyte isozyme) (Acylphosphatase;organ-common type isozyme) (Acylphosphate phosphohydrolase 1); erythrocyte isozyme) (Acylphosphatase; organ-common type isozyme) (Acylphosphate phosphohydrolase 1)
Observed Band	
Calculated Molecular Weight	11kD(Calculated).
Cell Pathway	
Tissue Specificity	Organ-common type isozyme is found in many different tissues.
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images