



SDF2L1 Mouse mAb

Catalog No	YP-mAb-18925
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SDF2L1 UNQ1941/PRO4424
Immunogen	Synthesized peptide derived from human SDF2L1
Specificity	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PWP1 interacting protein 8;SDF2 like protein 1;SDF2L1;SDF2L1 UNQ1941/PRO4424;Stromal cell-derived factor 2-like protein 1 (SDF2-like protein 1) (PWP1-interacting protein 8)
Observed Band	24 kDa
Calculated Molecular Weight	221 aa, 24 kDa
Cell Pathway	Endoplasmic reticulum lumen .
Tissue Specificity	Ubiquitously expressed with high expression in testis, moderate expression in the pancreas, spleen, prostate, small intestine and colon. Very low expression is seen in brain and skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images