



ZNF703 Mouse mAb

Catalog No	YP-mAb-18934
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF703 ZEPP01 ZPO1
Immunogen	Synthesized peptide derived from human ZNF703
Specificity	This antibody detects endogenous levels of ZNF703 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZPO1;ZNF503L;Zinc finger elbow-related proline domain protein 1;ZEPP01;ZNF703 ZEPP01 ZPO1;Zinc finger protein 703 (Zinc finger elbow-related proline domain protein 1)
Observed Band	58 kDa
Calculated Molecular Weight	590 aa, 58 kDa
Cell Pathway	Nucleus . Cytoplasm .
Tissue Specificity	Expressed in mammary epithelium.
Function	Transcriptional corepressor which does not bind directly to DNA and may regulate transcription through recruitment of histone deacetylases to gene promoters. Regulates cell adhesion, migration and proliferation. May be required for segmental gene expression during hindbrain development.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images