



TTC39A Mouse mAb

Catalog No	YP-mAb-18976
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TTC39A C1orf34 KIAA0452
Immunogen	Synthesized peptide derived from human TTC39A
Specificity	This antibody detects endogenous levels of TTC39A at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C1orf34;DEME 6;KIAA0452;TPR repeat protein 39A;TTC39A;TTC39A C1orf34 KIAA0452;Tetratricopeptide repeat protein 39A (TPR repeat protein 39A) (Differentially expressed in MCF7 with estradiol protein 6) (DEME-6)
Observed Band	65-70 kDa
Calculated Molecular Weight	613 aa, 70 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images