



PSG4 Mouse mAb

Catalog No	YP-mAb-18994
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PSG4 CGM4 PSG9
Immunogen	Synthesized peptide derived from human PSG4
Specificity	This antibody detects endogenous levels of PSG4 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CGM4;PSG9;PS beta G 4;PS beta G 9;PSBG 4;PSBG 9;PSG4;PSG4 CGM4 PSG9;Pregnancy-specific beta-1-glycoprotein 4 (PS-beta-G-4) (PSBG-4) (Pregnancy-specific glycoprotein 4) (Pregnancy-specific beta-1-glycoprotein 9) (PS-beta-G-9) (PSBG-9) (Pregnancy-specific glycoprotein 9)
Observed Band	
Calculated Molecular Weight	47 kDa
Cell Pathway	Secreted .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images