



KLHL35 Mouse mAb

Catalog No	YP-mAb-18995
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KLHL35
Immunogen	Synthesized peptide derived from human KLHL35
Specificity	This antibody detects endogenous levels of KLHL35 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	26597;kelch like 35 (Drosophila);Kelch like protein 35;KLHL35;Kelch-like protein 35
Observed Band	40-45 kDa
Calculated Molecular Weight	363 aa, 40 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images