



CORO2A Mouse mAb

Catalog No	YP-mAb-18996
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CORO2A IR10 WDR2
Immunogen	Synthesized peptide derived from human CORO2A
Specificity	This antibody detects endogenous levels of CORO2A at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CLIPINB;CORO2A;Coronin 2A;DKFZp686G19226;IR10;WD repeat containing protein 2;WDR2;CORO2A IR10 WDR2;Coronin-2A (IR10) (WD repeat-containing protein 2)
Observed Band	48 kDa
Calculated Molecular Weight	60 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images