



ATP6V0A4 Mouse mAb

Catalog No	YP-mAb-19014
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ATP6V0A4 ATP6N1B ATP6N2
Immunogen	Synthesized peptide derived from human ATP6V0A4
Specificity	This antibody detects endogenous levels of ATP6V0A4 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	A4;ATP6N1B;ATP6N2;ATP6V0A4;RDRTA2;RTA1C;RTADR;STV1;V ATPase 116 kDa isoform a4;VPH1;VPP2;ATP6V0A4 ATP6N1B ATP6N2;V-type proton ATPase 116 kDa subunit a isoform 4 (V-ATPase 116 kDa isoform a4) (Vacuolar proton translocating ATPase 116 kDa subunit a isoform 4) (Vacuolar proton translocating ATPase 116 kDa subunit a kidney isoform)
Observed Band	100 kDa
Calculated Molecular Weight	840 aa, 96 kDa
Cell Pathway	Apical cell membrane ; Multi-pass membrane protein . Basolateral cell membrane ; Multi-pass membrane protein . Localizes to the apical surface of alpha-intercalated cells in the cortical collecting ducts of the distal nephron (PubMed:10973252). Localizes to the basolateral surface of beta-intercalated cells in the cortical collecting ducts of the distal nephron (By similarity). .
Tissue Specificity	Expressed in adult and fetal kidney. Found in the inner ear.
Function	Subunit of the V0 complex of vacuolar(H ⁺)-ATPase (V-ATPase), a multisubunit enzyme composed of a peripheral complex (V1) that hydrolyzes ATP and a membrane integral complex (V0) that translocates protons (By similarity). V-ATPase is responsible for acidifying and maintaining the pH of intracellular compartments and in some cell types, is targeted to the plasma membrane,



where it is responsible for acidifying the extracellular environment (By similarity).
Involved in normal vectorial acid transport into the urine by the kidney .

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images