



STXBP6 Mouse mAb

Catalog No	YP-mAb-19054
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	STXBP6 HSPC156
Immunogen	Synthesized peptide derived from human STXBP6
Specificity	This antibody detects endogenous levels of STXBP6 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Amisyn;Syntaxin-binding protein 6;STXBP 6;HSPC156;STXBP6 HSPC156; Syntaxin-binding protein 6 (Amisyn)
Observed Band	28-30 kDa
Calculated Molecular Weight	25 kDa or 45 kDa
Cell Pathway	Cytoplasm . Membrane ; Peripheral membrane protein .
Tissue Specificity	Detected at low levels in brain, and at very low levels in heart, adrenal gland, testis, liver and kidney.
Function	Forms non-fusogenic complexes with SNAP25 and STX1A and may thereby modulate the formation of functional SNARE complexes and exocytosis.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images