



PLEKHA5 Mouse mAb

Catalog No	YP-mAb-19059
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PLEKHA5 KIAA1686 PEPP2
Immunogen	Synthesized peptide derived from human PLEKHA5
Specificity	This antibody detects endogenous levels of PLEKHA5 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA1686;PEPP 2;PEPP2;PLEKHA5;PLEKHA5 KIAA1686 PEPP2;Pleckstrin homology domain-containing family A member 5 (PH domain-containing family A member 5) (Phosphoinositol 3-phosphate-binding protein 2) (PEPP-2)
Observed Band	127 kDa, 74 kDa
Calculated Molecular Weight	1116 aa, 127 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	Highly expressed in heart and kidney.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images