



PRPSAP1 Mouse mAb

Catalog No	YP-mAb-19065
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PRPSAP1
Immunogen	Synthesized peptide derived from human PRPSAP1
Specificity	This antibody detects endogenous levels of PRPSAP1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	39 kDa phosphoribosypyrophosphate synthase-associated protein;PAP39; Phosphoribosyl pyrophosphate synthase-associated protein 1;PRPP synthase-associated protein 1;PRPSAP1;Phosphoribosyl pyrophosphate synthase-associated protein 1 (PRPP synthase-associated protein 1) (39 kDa phosphoribosypyrophosphate synthase-associated protein) (PAP39)
Observed Band	39 kDa
Calculated Molecular Weight	356 aa, 39 kDa
Cell Pathway	
Tissue Specificity	Ubiquitous.
Function	Seems to play a negative regulatory role in 5-phosphoribose 1-diphosphate synthesis.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images