



OTUD6A Mouse mAb

Catalog No	YP-mAb-19080
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	OTUD6A DUBA2
Immunogen	Synthesized peptide derived from human OTUD6A
Specificity	This antibody detects endogenous levels of OTUD6A at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DUBA 2;DUBA2;HSHIN6;OTU domain containing 6A;OTUD6A;OTUD6A DUBA2;OTU domain-containing protein 6A (DUBA-2)
Observed Band	33 kDa
Calculated Molecular Weight	288 aa, 33 kDa
Cell Pathway	
Tissue Specificity	
Function	Deubiquitinating enzyme that hydrolyzes 'Lys-27'-, 'Lys-29'- and 'Lys-33'-linked polyubiquitin chains. Also able to hydrolyze 'Lys-11'-linked ubiquitin chains.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images