



mon1b Mouse mAb

Catalog No	YP-mAb-19100
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MON1B HSRG1 KIAA0872 SAND2
Immunogen	Synthesized peptide derived from human MON1B
Specificity	This antibody detects endogenous levels of MON1B at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SAND2;KIAA0872;HSV-I stimulating-related protein;HSV-1 stimulation-related gene 1 protein;HSRG1;MON1B HSRG1 KIAA0872 SAND2;Vacuolar fusion protein MON1 homolog B (HSV-1 stimulation-related gene 1 protein) (HSV-I stimulating-related protein)
Observed Band	59 kDa
Calculated Molecular Weight	547 aa, 59 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images