



DTX4 Mouse mAb

Catalog No	YP-mAb-19110
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DTX4 KIAA0937 RNF155
Immunogen	Synthesized peptide derived from human DTX4
Specificity	This antibody detects endogenous levels of DTX4 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Deltex4;Dtx4;DTX4_HUMAN;E3 ubiquitin-protein ligase DTX4;Protein deltex-4; RING finger protein 155;RNF155;DTX4 KIAA0937 RNF155;E3 ubiquitin-protein ligase DTX4 (Protein deltex-4) (Deltex4) (RING finger protein 155)
Observed Band	
Calculated Molecular Weight	67 kDa; 67kD(Calculated).
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	Regulator of Notch signaling, a signaling pathway involved in cell-cell communications that regulates a broad spectrum of cell-fate determinations (By similarity). Functions as a ubiquitin ligase protein in vivo, mediating 'Lys48'-linked polyubiquitination and promoting degradation of TBK1, targeting to TBK1 requires interaction with NLRP4.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images