



PLAGL2 Mouse mAb

Catalog No	YP-mAb-19131
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PLAGL2 KIAA0198
Immunogen	Synthesized peptide derived from human PLAGL2
Specificity	This antibody detects endogenous levels of PLAGL2 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C2H2 type zinc finger protein;KIAA198;PLAGL2;PLAL2_HUMAN;Pleiomorphic adenoma gene-like 2;Pleiomorphic adenoma-like protein 2;Zinc finger protein PLAGL2;ZNF9;PLAGL2 KIAA0198;Zinc finger protein PLAGL2 (Pleiomorphic adenoma-like protein 2)
Observed Band	55kD
Calculated Molecular Weight	7 kDa; 55kD(Calculated).
Cell Pathway	Nucleus .
Tissue Specificity	
Function	Shows weak transcriptional activatory activity.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images