



CALML4 Mouse mAb

Catalog No	YP-mAb-19132
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CALML4
Immunogen	Synthesized peptide derived from human CALML4
Specificity	This antibody detects endogenous levels of CALML4 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CALML4;calmodulin like 4;Calmodulin like protein 4;NY BR 20;Calmodulin-like protein 4 (Serologically defined breast cancer antigen NY-BR-20)
Observed Band	22-25 kDa
Calculated Molecular Weight	196 aa, 22 kDa
Cell Pathway	
Tissue Specificity	Expressed in breast cancer cell lines.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images