



TTC17 Mouse mAb

Catalog No	YP-mAb-19150
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TTC17
Immunogen	Synthesized peptide derived from human TTC17
Specificity	This antibody detects endogenous levels of TTC17 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TPR repeat protein 17;TTC17;Tetratricopeptide repeat protein 17 (TPR repeat protein 17)
Observed Band	
Calculated Molecular Weight	1141 aa, 13 kDa
Cell Pathway	Cytoplasm . Cell membrane . Cytoplasm, cytoskeleton . Colocalized with CATIP at F-actin rich zones and at dynamic plasma membrane protrusions.
Tissue Specificity	Expressed in germ cells as well as in somatic cells of the testis (at protein level).
Function	Plays a role in primary ciliogenesis by modulating actin polymerization.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images